

HOW TO ATTACH YOUR



HAPPY BACK GEL CELL LUMBAR CUSHION

Item # 2606-R



happy BACK



Follow these SIMPLE STEPS

CONTACT US!  
We're here to help

✉ info@novajoy.com

☎ 1 800 557 6682

novajoy.com

Method 1:

To attach the straps (crisscross)

- 1 Position the cushion on the back of the chair so the slip resistant back is touching the chairs back. Un-buckle the straps and place around back of chair. The two top straps go over the back and the two bottom go around the back (see Figure 2).



Figure 1

- 2 Adjust the length of the strap by sliding the plastic slider to shorten or lengthen the strap (see Figure 2).



PLASTIC SLIDER Figure 2

- 3 Buckle the straps together making a crisscross as shown below blue and green arrows (see Figure 3).  
NOTE: Ensure the straps are not loose.



Figure 3

Method 2:

To attach the straps (horizontally)

- 1 Position the cushion on the back of the chair so the slip resistant back is touching the chairs back. Un-buckle the straps and place around back of chair (see Figure 1).




Figure 1

- 2 Adjust the length of the strap by sliding the plastic slider to shorten or lengthen the strap (see Figure 2).



PLASTIC SLIDER Figure 2

- 3 Buckle the straps together as shown below blue and green lines (see Figure 3).  
NOTE: Ensure the straps are not loose.



Figure 3

Method 3:

To attach the straps (up & around)

- 1 Position the cushion on the back of the chair so the slip resistant back is touching the chairs back. Un-buckle the straps and place around back of chair (see Figure 1).



Figure 1

- 2 Adjust the length of the strap by sliding the plastic slider to shorten or lengthen the strap (see Figure 2).



PLASTIC SLIDER Figure 2

- 3 Buckle the straps together as shown below blue and green lines (see Figure 3).  
NOTE: Ensure the straps are not loose.



Figure 3

