



nova

353

Parts Breakdown
For models: 353B, 353BK, 353P, and 353R



An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your transport chair corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (see Figure 1).

If you need further assistance, please call our customer service team at (800) 557-6682.

Thank you!



Figure 1

Hand Brake Mechanism

CH-35307



CH-33016



CH-35304



Handle and Folding Mechanism Components

CH-35313

CH-35315

CH-35309

CH-35311

CH-35312

CH-35314

CH-35308

CH-35310



CH-37701



CH-37724



CHT-31954



CH-37715



CH-35003

CH-35008

CH-35002

CH-35007

CH-35303

CH-35317

CH-35302

CH-35316



CH-35018



CH-35005




Main Body Components




Footrest Assembly



Front Wheel Assembly



Rear Wheel Assembly



	Item Number 353
	353 Transport Chair
	SERIAL # BEGINNING WITH "CH"
ITEM NUMBER	DESCRIPTION
CH-35301	ANTI-TIPPER (PAIR) BLACK
CH-37701	ARM PAD FOR 377/379 SN:CH
CH-33014	AXLE & NUT FOR REAR WHEEL
CH-31925	AXLE SET FOR FRONT WHEEL 309
CH-33006	BEARING & SPACER (SET) FOR REAR
CH-31929	BEARING & SPACER FOR FRONT WHEEL
CH-31933	BEARINGS FOR FRONT FORK 307
CH-35302	BLACK LEFT D/A FOR 353
CH-35303	BLACK RIGHT D/A FOR 353
CH-35002	BLUE LEFT D/A FOR 350
CH-35003	BLUE RIGHT D/A FOR 350
CH-33012	BRAKE LEFT FOR 330 SN:CH
CH-33013	BRAKE RIGHT FOR 330 SN:CH
CH-37728	CAP FOR FOOT REST SN: CH
CH-35304	COVER FOR BRAKE CABLE FOR 353
PVT-AC101	COVER FOR REAR WHEEL AXLE (330)
CH-35004	CRADLE FOR FRONT SEAT
CH-35005	DESK ARM PANEL W/ HARDWARE LEFT
CH-35018	DESK ARM PANEL W/ HARDWARE RIGHT
CH-30707 (PAIR)	FOOT PLATE (PAIR) FOR 307/309
CH-35305	FOOTRESTS BLACK (PAIR) FOR 327/329/330
CHT-32710 (PAIR)	FOOTRESTS BLUE (PAIR) FOR 327/329/330
CH-35306	FOOTRESTS PURPLE (PAIR) FOR 327/329/330
CHT-32711 (PAIR)	FOOTRESTS RED (PAIR) FOR 327/329/330
CH-31906	FORK FOR FRONT WHEEL 319/327/330
CH-33016 (PAIR)	GRIP (PAIR) FOR 330 SN: CH
CH-35307	HAND BRAKE W/CABLE FOR 353
CH-35308	HANDLE LEFT FOR 353 BLACK
CH-35309	HANDLE LEFT FOR 353 BLUE
CH-35310	HANDLE LEFT FOR 353 PURPLE
CH-35311	HANDLE LEFT FOR 353 RED
CH-35312	HANDLE RIGHT FOR 353 BLACK
CH-35313	HANDLE RIGHT FOR 353 BLUE
CH-35314	HANDLE RIGHT FOR 353 PURPLE
CH-35315	HANDLE RIGHT FOR 353 RED
CH-30714	HEEL STRAP SN: CH
CH-37718	KNOB WITH SCREW FOR FOOTREST

	Item Number 353
	353 Transport Chair
	SERIAL # BEGINNING WITH "CH"
ITEM NUMBER	DESCRIPTION
CH-37724	LEVER FOR FOLDING MECHANISM
CH-37717 (PAIR)	PIN FOR FOOTREST
CH-31934	PLUG FOR NUT FOR FRONT FORK
CH-35316	PURPLE LEFT D/A FOR 353
CH-35317	PURPLE RIGHT D/A FOR 353
CH-35006	RECEIVER FOR FLIP BACK ARM-REAR FOR 350 W/ HARDWARE
CH-35007	RED LEFT D/A FOR 350
CH-35008	RED RIGHT D/A FOR 350
CH-35009	RELEASE LEVER ON FLIP BACK ARM
CH-37737	SCREW FOR BRAKE SHOE
CH-30713	SCREW FOR HEEL STRAP SN: CH
CH-37730	SCREW FOR UPHOLSTERY/HAND BRAKE
CT-SB100	SEAT BELT FOR ALL TRANSPORT CHAIRS
CH-31935	TIP FOR BRAKE HANDLE 307/309/3
CH-35318	UPHOLSTERY BACK FOR 353
CH-31903	UPHOLSTERY SEAT FOR 330 SN:CH
CH-37731	WASHER FOR UPHOLSTERY SCREW
CH-33003	WHEEL FRONT FOR 330 SN:CH
CH-33004	WHEEL REAR FOR 330 SN:CH

CONTACT US!
We're here to help



info@novajoy.com



1 800 557 6682

novajoy.com

CARE • QUALITY • INNOVATION



THE **NOVA** PROMISE