



Mini Mack

Parts Breakdown
For models: 4214BL, 4214RD
Serial # Starting "NT"



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

Revised on February 5, 2019

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (**see Figure 1**).


If you need further assistance, please call our customer service team at **(800) 557-6682**.

Thank you!




Figure 1


Hand Brake Mechanism



Upper Assembly Components




Seat & Backrest



Folding Mechanism & Accessories



Wheel & Lower Assembly Components






Item Number 4214

Mini Mack


SERIAL # BEGINNING WITH "NT"

ITEM NUMBER	DESCRIPTION	UOM
P42004	ANTI-RATTLE SEAT CLAMP FOR PLASTIC SEAT	PAIR
P40021	AXLE FEMALE FOR REAR WHEEL	EACH
P42090	AXLE FOR FRONT WHEEL	EACH
P40017	AXLE REAR WHEEL - MALE	EACH
P421515BL	BACKREST BLUE FLIP-UP	EACH
P421515RD	BACKREST RED FLIP-UP	EACH
P4214B1	BASKET	EACH
VT-33207	BEARING FOR SWIVEL FRONT FORK	EACH
P40006	BEARINGS FOR WHEEL	BAG5
P-CABLE	BRAKE CABLE (PAIR)	PAIR
P40027	BRAKE SHOE	PAIR
P-COVER	CABLE COVER BLACK	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
V42093	FOAM FOR FLIP-UP BACKREST	EACH
P42191FM	FOLDING MECHANISM	EACH
KP4218F	FORK ASSEMBLY FOR 8" FRONT WHEEL	EACH
P40033	FORK FOR 8" FRONT WHEEL W/STEM	EACH
P40031	FORK FOR REAR WHEEL	EACH
P43003	GRIP - LEFT ANATOMICAL	EACH
P43004	GRIP - RIGHT ANATOMICAL	EACH
B14-C	HAND BRAKE WITH STEEL CABLE	EACH
P42159	HANDLE - LEFT	EACH
P421510	HANDLE - RIGHT SIDE	EACH
P421512	HANDLE SILENCER	PAIR
V420902	HARDWARE FOR FOLDING MECHANISM	EACH
H-1002	LOCK N LOAD	EACH
P40015	NUT COVER FOR WHEEL AXLES	EACH
P40028	NUT FOR FOLDING MECH	EACH
P40011	NUT FOR FORK STEM	EACH
P40016	NUT FOR FRONT/REAR WHEEL AXLE	BAG5
P42058	PAD FOR SEAT	EACH
P-421401	PLUG FOR SLANTED FOLDING FRAME	EACH
V42076	RECEIVER FOR REMOVABLE BACKREST	EACH
P40047	REFLECTOR FOR BRAKE HANDLE	EACH
V42064	SCREW & NUT FOR BACKREST HINGE	EACH
P40045	SCREW FOR FOLDING MECHANISM	EACH

	Item Number 4214	
	Mini Mack	
	SERIAL # BEGINNING WITH "NT"	
ITEM NUMBER	DESCRIPTION	UOM
P42007	SCREW FOR FOLDING MECH. & SLANTED FRAME PLUG	EACH
V42009	SCREW SET FOR A-FRAME	EACH
P40007	SCREW TO ADJUST CABLE TENSION	EACH
P40030	SCREW TO ATTACH BRAKE SHOE TO REAR FORK	EACH
P40035	SCREW TO ATTACH FORK TO FRAME	EACH
C40006	SCREW TO ATTACH HAND BRAKE	EACH
P40029	SCREW TO ATTACH REAR FORK TO FRAME	EACH
P40001	SCREWS TO ATTACH PLASTIC SEAT	BAG5
P4215S	SEAT	EACH
V42089	TRIANGLE SHAPED HEIGHT ADJUSTMENT	EACH
P42003	U-SHAPE SEAT ATTACHMENT	EACH
P42009	WASHER FOR FOLDING MECHANISM CONTOURED	EACH
P42010	WASHER PLASTIC CONTOURED FOR FOLDING MECHANISM	EACH
P421517	WASHERS FOR WHEEL	PAIR
P421516	WHEEL	EACH
P42003	U-SHAPE SEAT ATTACHMENT	EACH
P42009	WASHER FOR FOLDING MECHANISM	EACH
P42010	WASHER PLASTIC CONTOURED FOR FOLDING MECH	EACH
P42085	WHEEL 8"	EACH

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

CARE • QUALITY • INNOVATION



THE NOVA PROMISE