



Forté 23

Parts Breakdown
*For models: 4323BL, 4323DB, 4323RD
Serial # Starting "NV"*



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

Revised on February 5, 2019

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (**see Figure 1**).


If you need further assistance, please call our customer service team at **(800) 557-6682**.

Thank you!



Figure 1

Hand Brake Mechanism




Upper Assembly Components




Seat & Backrest



Folding Mechanism And Toe Kick



Wheel & Lower Assembly Components






Item Number 4323

Forté 23


SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
V-432002	AXLE FRONT WITH NUT	SET
V-432003	AXLE REAR FEMALE AND MALE SCREW	SET
V-432014	BACK REST	EACH
V-432026	BAG	EACH
V-432032	BAR (2PCS) FOR SEAT	PAIR
V-432004	BEARINGS FOR 8" WHEEL WITH SPACERS (PAIR)	PAIR
V-432031	BEARINGS WITH SPACER FOR FRONT SWIVEL FORK	SET
V-432001	BRAKE ASSEMBLY W/ CABLE & COVER	EACH
P-CABLE	BRAKE CABLE (PAIR)	PAIR
P-TKC015	BRAKE CABLE HOLDER	EACH
V42087	BRAKE SHOE	EACH
JN-490004	COVER FOR CABLE	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
V-432036	CROSSBAR SUPPORT (PAIR)	PAIR
V-432035	CROSSBAR SUPPORT (PAIR)	EACH
V-432019	FORK ATTACHMENT WITH SCREW	EACH
V-432015	FORK BOTTOM HALF WITH SWIVEL BEARING	EACH
V-432033	FRAME ATTACHMENT LEFT	EACH
V-432034	FRAME ATTACHMENT RIGHT	EACH
V-432010	GRIPS (PAIR)	PAIR
V-432011	HANDLE BAR FOR THE LEFT SIDE	EACH
V-432012	HANDLE BAR FOR THE RIGHT SIDE	EACH
V-432017	KNOB LOCKING RELEASE FOR FOLDING MECH	EACH
P-TKW1212	KNOB WITH SCREW FOR HANDLE HEIGHT	EACH
P42012	NUT FOR BRAKE SHOE SCREW	EACH
V-432008	PLUG FOR SEAT (PAIR)	PAIR

	Item Number 4323	
	Forté 23	
	SERIAL # BEGINNING WITH "NV"	
ITEM NUMBER	DESCRIPTION	UOM
P40047	REFLECTOR FOR BRAKE HANDLE	PAIR
P40103	SCREW & NUT (PAIR) TO ATTACH BRAKE HANDLE	PAIR
P42015	SCREW FOR BRAKE SHOE	EACH
V-432028	SCREW SET FOR BACKREST	PAIR
V-432021	SCREW SET FOR CENTER OF FOLDING MECH	SET
V-432020	SCREW SET FOR CROSSBAR SUPPORT	PAIR
V-432016	SCREW SET FOR FRONT FORK	SET
P40007	SCREW TO ADJUST CABLE TENSION	EACH
V-432009	SEAT	EACH
P42013	SPRING FOR BRAKE (PAIR)	PAIR
V-432013	TOE KICK (PAIR)	PAIR
V42074D	WHEEL 8"	EACH
P40014	WHITE PLASTIC WASHER FOR BRAKE	EACH

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

CARE • QUALITY • INNOVATION



THE **NOVA** PROMISE