



4500 gemini

Two Rides ... One Destination

Parts Breakdown
For model: 4500BK



An important note from




Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.


Please verify that the serial number on your unit corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (see Figure 1).

If you need further assistance, please call our customer service team at (800) 557-6682.

Thank you!



Hand Brake Mechanism



Main Body Components



Footrest Assembly



Front Wheel Assembly

RM-450015



RM-450016



RM-450017



RM-450025



Rear Wheel Assembly

RM-450018



RM-450027



RM-450019



RM-450008



RM-450007



	Item Number 4500
	4500 GEMINI 2 in 1 Rollator Transport Chair Combo
	SERIAL # BEGINNING WITH "RM"
ITEM NUMBER	DESCRIPTION
RM-450001	POUCH FOR 4500
RM-450002	BACKREST FOR 4500 W/HARDWARE
RM-450003	SEAT FOR 4500
RM-450004	FOOTRESTS (PAIR) FOR 4500
RM-450005	HANDLE ASSEMBLY LEFT
RM-450006	HANDLE ASSEMBLY RIGHT
RM-450007	BRAKE SHOE LEFT W/HARDWARE
RM-450008	BRAKE SHOE RIGHT W/HARDWARE
RM-450009	SEAT BELT FOR 4500
RM-450010	FOOT PLATE FOR LEFT FOOTREST
RM-450011	HEEL STRAP for 4500
RM-450012	ARM PAD LEFT W/HARDWARE
RM-450013	ARM PAD RIGHT W/ HARDWARE
RM-450014	GRIP for 4500
RM-450015	FRONT WHEEL W/ BEARING AND SPACER
RM-450016	BEARING AND SPACER FOR FRONT WHEEL
RM-450017	FRONT FORK ASSEMBLY FOR 4500
RM-450018	REAR WHEEL W/ BEARING AND SPACER
RM-450019	AXLE FOR REAR WHEEL
RM-450020	CABLE COVER
RM-450021	FRONT FORK BEARINGS
RM-450022	HANDLE HEIGHT ADJUSTMENT KNOB
RM-450023	SEAT SILENCER
RM-450024	SEAT HANDLE WITH KEY RING
RM-450025	AXLE FOR FRONT WHEEL
RM-450026	FOOT PLATE FOR RIGHT FOOTREST
RM-450027	BEARING AND SPACER FOR REAR WHEEL