



# nova

## 5240

Parts Breakdown  
*For models: 5240S, 5240SE, 5241S, 5241SE  
Serial # Starting "YH"*



CARE • QUALITY • INNOVATION



THE NOVA PROMISE

*Revised on December 10, 2018*

**CONTACT US!**  
We're here to help

 [info@novajoy.com](mailto:info@novajoy.com)

 1 800 557 6682

[novajoy.com](http://novajoy.com)

## An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your wheelchair corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame near the front wheel (**see Figure 1**).


If you need further assistance, please call our customer service team at **(800) 557-6682**.

**Thank you!**




**Figure 1**


## Main Body Components



## Footrest & Lower Frame Components





## Front Wheel Assembly






## Rear Wheel Assembly




	Item Number 5240S/5240SE/5241S/5241SE	
	5240 Steel Wheelchair	
	SERIAL # BEGINNING WITH "YH"	
ITEM NUMBER	DESCRIPTION	UOM
YH-616018	ANTI-TIPPER (PAIR)	PAIR
YU-506019	ARM PAD DESK ARM WITH SCREWS	EACH
YU-506052	ARM PAD FULL LENGTH WITH SCREW	EACH
YU-522004	AXLE REAR	EACH
YU-506015	AXLE SET FOR 8" WHEEL	SET
YU-522003	BEARING & SPACER SET FOR 24" W	SET
YU-506012	BEARING & SPACER SET FOR 8" WHEEL	SET
YH-522001	BRAKE LEFT	EACH
YH-522002	BRAKE RIGHT	EACH
YU-506011	CAP FOR FOOTREST (PAIR)	PAIR
YU-522007	CRADLE FOR 5220S/5240S	EACH

	Item Number 5240S/5240SE/5241S/5241SE	
	<b>5240 Steel Wheelchair</b>	
	SERIAL # BEGINNING WITH "YH"	
ITEM NUMBER	DESCRIPTION	UOM
YH-5100LD1	DESK ARM LEFT	EACH
YH-5100RD1	DESK ARM RIGHT	EACH
W-7001	ELEVATING LEG RESTS	PAIR
YU-522006	FOOT PLATE METAL (PAIR)	PAIR
YH-520002	FOOTRESTS BLACK	PAIR
YU-506045	FORK FOR FRONT WHEEL	EACH
YH-5101LD1	FULL ARM LEFT	EACH
YH-5101RD1	FULL ARM RIGHT	EACH
YU-506039	GRIP (PAIR)	PAIR
YU-522008	HAND RIM CHROME FOR 24" WHEEL	EACH
YU-506053	HUB CAP FOR 24" WHEEL	PAIR
YU-516003	PANEL LEFT FOR DESK	EACH
YU-516004	PANEL RIGHT FOR DESK	EACH
YH-616007	PANEL LEFT FOR FULL ARM	EACH
YH-616008	PANEL RIGHT FOR FULL ARM	EACH
YU-506054	PLUG (PAIR) FOR FOOT REST	PAIR
YU-506048	PLUG (PAIR) FOR FRAME	PAIR
YU-506024	PIN (4PCS) FOR FOOTREST	SET
YU-522009	SCREW CENTER FOR CROSSBAR	EACH
YU-506055	SCREW FOR ARM PAD	PAIR
YU-506058	SCREW FOR PANEL/SKIRT GUARD	PAIR
YU-506034	SEAT BELT	EACH
S5100SB	SWIVEL BEARING (PAIR)	PAIR
YU-506051	TIP FOR BRAKE	EACH
YU-524001	UPHOLSTERY BACK	EACH
YU-524002	UPHOLSTERY SEAT	EACH
S5100FW2	WHEEL 8" FRONT	EACH
YU-522005	WHEEL 24" REAR	EACH

**CONTACT US!**  
We're here to help

 [info@novajoy.com](mailto:info@novajoy.com)

 1 800 557 6682

[novajoy.com](http://novajoy.com)

CARE • QUALITY • INNOVATION



THE NOVA PROMISE