



TKW-14

Parts Breakdown
*For models: TKW-14
Serial # Starting "NT"*



CARE • QUALITY • INNOVATION



THE NOVA PROMISE

Revised on October 15, 2024

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

[novajoy.com](https://www.novajoy.com)

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your Knee Walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker on the side of the frame (see Figure 1).


If you need further assistance, please call our customer service team at (800) 557-6682.

Thank you!




Figure 1



Steering And Brake Components



Front And Rear Wheel Components



Knee Pad Platform And Accessories





Item Number TKW-14

THE KNEE WALKER

SERIAL # BEGINNING WITH "NT"

ITEM NUMBER	DESCRIPTION	UOM
P-TKW1411	AXLE FOR WHEEL	EACH
P-TKW1215	BARREL SCREW SET (PAIR) FOR BRAKE	PAIR
TKW-BASKET	BASKET	EACH
P-TKW1412	BEARING FOR WHEEL	SET
P-TKW1405	BEARING FOR STEERING	PAIR
P-TKW1207	BOLT SET FOR FOLDING HINGE	EACH
P-TKW1201	BRAKE LEFT	EACH
P-TKW1227	BRAKE PAD (PAIR)	PAIR
P-TKW1202	BRAKE RIGHT	EACH
P-TKW1214	CABLE W/COVER FOR BRAKE	EACH
P-TKW1406	C-CLAMP FOR STEERING	EACH
P40013	CRIMP FOR STEEL CABLE	PAIR
P-TKW1229	GRIP - PAIR	PAIR
B15-C	HAND BRAKE WITH STEEL CABLE SEE SN REFERENCE BELOW	EACH
P-TKW1213	HAND BRAKE W/CABLE	EACH
V42010	HAND GRIPS ANATOMICAL SEE SN REFERENCE BELOW	PAIR
P-TKW1415	HANDLE BAR SEE SN REFERENCE BELOW	EACH
P-TKW1410	HANDLE BAR SEE SN REFERENCE BELOW	EACH
P-TKW1414	HANDLE BAR	EACH
P-TKWBH	HANDLE FOR TKW-BASKET	EACH
P-TKW1224	HEIGHT ADJUSTMENT LEVER (ORANGE)	SET
P-TKW1221	HOOK FOR BASKET	SET
P-TKW1209	KNEE PAD	EACH
P-TKW1404	KNEE PLATFORM	EACH
P-TKW1212	KNOB WITH SCREW FOR HANDLE HEIGHT SEE SN REFERENCE BELOW	EACH
P-TKW1206	LEVER FOR FOLDING MECHANISM	SET
P-TKW1408	PLUG REAR FRAME	EACH
P-TKW1402	POST LOWER STEERING	EACH
P40103	SCREW & NUT TO ATTACH BRAKE HANDLE	PAIR
P-TKW1203	SCREW & WASHER FOR BRAKE	SET
P-TKW1217	SCREW FOR BRAKE ADJUSTMENT (PAIR)	PAIR
P-TKW1407	SCREW FOR LOWER STEERING POST	EACH

	Item Number TKW-14	
	THE KNEE WALKER	
	SERIAL # BEGINNING WITH "NT"	
ITEM NUMBER	DESCRIPTION	UOM
P-TKW1211	SCREW SET (8 PCS) FOR SEAT	SET
P-TKW1208	SCREW W/ KNOB FOR HANDLE BAR AND SEAT HT. ADJ	EACH
P-TKW1409	SILENCER FOR HANDLE	EACH
P-TKW1216	SPRING (PAIR) FOR BRAKE	PAIR
P-TKW1403	STEERING POST SEE SN REFERENCE BELOW	SET
P-TKW1413	STEERING POST	EACH
P-TKW1401	WHEEL	EACH
P-TKW1220	WRENCH (2 PCS) FOR TKW-12	EACH

SERIAL NUMBER (SN) REFERENCE	
ITEM NUMBER	<i>PRE JAN. 2017 MODELS</i>
TKW-14	SN BEFORE: NTTKW-14FN0001

SERIAL NUMBER (SN) REFERENCE	
ITEM NUMBER	<i>POST MAR. 2020 MODELS</i>
TKW-14	SN BEFORE: NTTKW-14KB0001