

## Limited Lifetime Warranty

Nova stands firm on our commitment to providing the highest quality products, as well as exceptional service. This product was built to precise standards and thoroughly inspected prior to shipment. This warranty represents our confidence in the materials and workmanship of our product.

Nova extends this warranty only to the original purchaser of this product. The warranty does not extend to any subsequent purchaser or owner. The warranty is void upon any sale or transfer of ownership or use by other person. Nova warrants its products are free from defects in material and workmanship under NORMAL use for the period commencing upon the date of purchase and continuing for the following specified period of time after that date:

Metal frame: Limited Lifetime Warranty

Non-Metal/Metal parts: One year Limited Warranty  
(wheels, hand grips, rubber tips, plastic)

This warranty does not cover problems due to user negligence, misuse of the product, failure to adhere to product instructions, or normal wear and tear. Items that have been subject to negligent abuse, modified without prior written authorization are not covered within this warranty.

For warranty service, please contact the dealer from whom you purchased your Nova product. In the event that you do not receive satisfactory service, please call us directly at the number below. Be prepared to indicate, the nature of the defect, the product serial number, and the name and location of where you purchased the product. You will also be required to provide a receipt of your purchase. DO NOT return any products to our office without prior consent.

Repair or replacement is the only remedy under this limited warranty. This warranty does not include any labor or shipping charges incurred in replacement part installation or repair of product.

THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER EXPRESS WARRANTIES, IMPLIED WARRANTIES, IF ANY, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NOVA SHALL NOT BE LIABLE FOR ANY SEQUENTIAL OR INCIDENTAL DAMAGES WHATSOEVER.

Customer help line:

**1-800-557-6682**

www.novamedicalproducts.com  
P.O. Box 3039 • Gardena CA 90247



## Travel Seat Cane with Sling Seat

---

Easy to use

Great for traveling

---

Products for an Independent Life.

www.novamedicalproducts.com



# How to use your NOVA Travel Seat Cane with Sling Seat

Item # 3092

Capacity not to exceed 250 lbs.

## For Walking:

- 1 Fold seat cane.
- 2 Wrap seat upholstery around both middle crossbars and secure with velcro strap.
- 3 Use the cane in the hand of the opposite side of your injured or weaker leg.

## For Sitting:

- 1 Unfold seat cane so the sling seat is completely flat. (See image 1)
- 2 Place on a flat level surface.
- 3 Carefully sit down on the seat with the larger handle in front. (See image 2)
- 4 Hold on to handle in front for additional stability. (See image 2)

Image 1





Image 2



Use only on a flat level surface.

Seat must be flat with all 4 legs touching the ground.

Always sit facing the larger handle and hold handle when seated.

Serious injury may occur if seating instructions are not followed.