

PLEASE READ BEFORE SEARCHING FOR THE PART NEEDED



Many of the parts pictured on this breakdown may look identical to the part or parts needed but the measurements may be different and will not be interchangeable. Please check the serial number on your transport chair and make sure it corresponds with the serial number on this part diagram highlighted in yellow. The serial numbers on our transport chairs are located on a sticker on the crossbar under the seat as shown in the picture below. If you cannot find the serial number, please call customer service for assistance at (800) 557-6682 option 1.

Thank you



NOVA Customer Service: 800-557-6682

www.NovaMedicalProducts.com

ITEM NUMBER 377 LIGHTWEIGHT ALUMINUM TRANSPORT CHAIR

SERIAL NUMBER BEGINNING WITH "CH"



CH-37706



CH-30710



CH-31935



CH-37701



CH-37702



CH-37703



CH-37704



CH-37705



CH-37708



CH-37709



CH-37710



CH-37711



CH-37712



CH-37713



CH-37714



CH-37715



CH-37716



CH-37717



CH-37718




CH-37719



CH-37721



CH-37723



CH-37724



CH-37728 CH-37729 CH-37730 CH-37731 CH-37732



CH-37733



CH-37734



CH-37735



CH-37736 CH-37737 CH-37738



CH-37739



CH-37740 CH-37741



CH-T1000




CT-SB100



NOVA Customer Service: 800-557-6682

www.NovaMedicalProducts.com

ITEM NUMBER 377	LIGHTWEIGHT ALUMINUM TRANSPORT CHAIR	
	SERIAL NUMBER BEGINNING WITH "CH"	
ITEM/PART NUMBER	DESCRIPTION	
CH-37706	WHEEL FRONT	
CH-30710	WHEEL REAR	
CH-31935	TIP FOR BRAKE HANDLE	
CH-37701	ARM PAD (INCLUDES SCREWS)	
CH-37702	SEAT UPHOLSTERY BLACK	
CH-37703	BACK UPHOLSTERY BLACK	
CH-37704	STAR KNOB FOR FRONT FORK AND REAR LEG	
CH-37705	FORK FRONT (OLD STYLE WITH SILVER NUT ON STEM)	
CH-37708	BRAKE LEFT (INCLUDES SCREW, WASHER & NUT)	
CH-37709	BRAKE RIGHT (INCLUDES SCREW, WASHER & NUT)	
CH-37710	FOLDING MECHANISM	
CH-37711	HANDLE LEFT BLUE	
CH-37712	HANDLE RIGHT BLUE	
CH-37713	GRIP	
CH-37714	TOE KICK	
CH-37715	SCREW FOR HANDLE & FOLDING MECHANISM	PAIR
CH-37716	FOOT PLATE	
CH-37717	PIN FOR FOOT REST	PAIR
CH-37718	KNOB FOR FOOT REST	
CH-37719	FOOT REST BLUE	PAIR
CH-37721	AXLE FOR FRONT WHEEL (INCLUDES NUT)	PAIR
CH-37723	BEARING & SPACER FOR FRONT WHEEL	PAIR
CH-37724	LEVER FOR FOLDING MECHANISM	
CH-37728	CAP FOR FOOT REST	
CH-37729	RELEASE PIN FOR FOLDING MECHANISM	
CH-37730	SCREW FOR UPHOLSTERY	
CH-37731	WASHER FOR UPHOLSTERY	SET
CH-37732	AXLE FOR REAR WHEEL (INCLUDES WASHERS & NUT)	PAIR
CH-37733	BEARING & SPACER FOR REAR WHEEL	PAIR
CH-37734	CROSSBAR SUPPORT	
CH-37735	SCREWS FOR CROSSBAR SUPPORT (INCLUDES WASHERS & NUTS)	PAIR
CH-37736	SCREW FOR CROSSBAR (INCLUDES WASHERS, NUT & SPACERS)	
CH-37737	SCREW, WASHER, AND NUT FOR BRAKE	SET
CH-37738	FORK FRONT NEW STYLE WITHOUT SILVER NUT ON STEM)	
CH-37739	PLUG, SPING, AND PIN FOR FRONT FORK AND REAR LEG	
CH-37740	REAR LEG ATTACHMENT LEFT	
CH-37741	REAR LEG ATTACHMENT RIGHT	

CH-T1000	TOOLS	
CT-SB100	SEAT BELT	

NOVA Customer Service: 800-557-6682

www.NovaMedicalProducts.com