

CRUISER II- ITEM # 4201



40023BK



40026GR C4203B



KP4200F



P40003 P40006 P40007 P40016 P40028



V42065



P40032



P42004 P42007 P42009 P42010 P42019 P42022



P42046 P43003 P43004 V42005 V42006 V42007 V42031 V42064



V42073




V42076 V42089 V42093



V42094



ITEM NUMBER 4201	CRUISER II	
	Example SN# 1104010145V4201AA0000019	
ITEM NUMBER	DESCRIPTION	
40023BK (pair)	TIPS FOR 7/8in. WALKER (4200,4200C,4201)	\$9.95
40026GR	WALKER GLIDE SKIS FOR 1 in. WALKER	\$19.95
C4203B	BASKET FOR 4200/4200C/4202/4202C/4203/4207/4212	\$19.95
KP4200F	FORK ASSEMBLY FOR 5" FRONT WHEEL (P40016, P42022, P42046)	\$9.95
P40003	SILENCER FOR ADJ. HANDLE	\$2.50
P40006 (pack 5)	BEARINGS FOR 5"/6"/8" WHEELS	\$9.95
P40007	SCREW TO ADJUST CABLE TENSION	\$2.50
P40016	NUT FOR FRONT/REAR WHEEL AXLE	\$2.50
P40028	STEEL NUT FOR FOLDING MECH	\$2.50
P40032	FOLDING MECHANISM FOR V42073	\$6.95
P42004 (pair)	ANTI-RATTLE SEAT CLAMP FOR PLASTIC SEAT	\$2.50
P42007	SCREW FOR FOLDING MECH. & SLANTED FRAME PLUG	\$3.50
P42009	STEEL CONTOURED WASHER FOR FOLDING MECHANISM	\$2.50
P42010	PLASTIC CONTOURED WASHER FOR FOLDING MECHANISM	\$2.50
P42019	WHEEL 5" FOR 4200/4201/4010 (INCLUDES BEARINGS)	\$9.95
P42022	FORK FOR 5" FRONT WHEEL	\$9.95
P42046	AXLE FOR FRONT WHEEL	\$2.50
P43003	LEFT BLACK ANATOMICAL GRIP	\$3.50
P43004	RIGHT BLACK ANATOMICAL GRIP	\$3.50
V42005	LEFT ALUM. HANDLE W/ GRIP FOR 4200/4200C	\$9.95
V42006	RIGHT ALUM. HANDLE W/ GRIP FOR 4200/4200C	\$9.95
V42007	SEAT CLAMP Fits 4200/4202C/4203/4208/4212/ALL ZOOMS	\$2.50
V42031	COMPLETE HARDWARE FOR V42073	\$3.50
V42064	SCREW & NUT FOR HINGE FOR BACKREST	\$2.50
V42065	BLUE DETACHABLE FLIP UP BACKREST	\$19.95
V42073	SEAT FOR - (4212/4208/4203/4201/4200)	\$19.95
V42076	BLACK PLASTIC HINGE FOR BACKREST	\$2.50
V42089	TRIANGLE SHAPED HEIGHT ADJUSTMENT KNOB	\$6.95
V42093	FOAM CUSHION FOR FLIP UP BACKREST	\$6.95
V42094	SLANTED FRAME FOLDING PLUG	\$2.50

NOVA Customer Service: 800-557-6682  
[www.NovaMedicalProducts.com](http://www.NovaMedicalProducts.com)