



GETGO- ITEM # 4203



V42064 V42073



V42073H



V42076 V42079 V42080



V42089 V42093



V42094 V42097



EXAMPLE SERIAL NUMBER 1106040203V4206AA0004551
AND SERIAL NUMBERS BEGINNING WITH "NV"

NOVA Customer Service: 800-557-6682
www.NovaMedicalProducts.com

ITEM NUMBER 4203		
	GETGO	
	SERIAL NUMBERS BEGINNING WITH "NV"	
ITEM NUMBER	DESCRIPTION	MSRP
B14-C	B14 HAND BRAKE WITH STEEL CABLE	\$19.95
C40006	SCREW TO ATTACH B11/B14 TO HANDLE	\$3.50
C4203B	WIRE BASKET	\$19.95
H-1002	LOCK N LOAD	\$19.95
KP4202F	FORK ASSEMBLY FOR 6" FRONT WHEEL	\$19.95
	(INCLUDES P42021, P40030, P40016)	
P40001 (5 PACK)	SCREWS FOR U-SHAPED SEAT ATTACHMENTS AND SEAT	\$2.50
P40003	HANDLE SILENCER	\$2.50
P40006 (5 PACK)	BEARINGS FOR 5"/6"/8" WHEELS	\$9.95
P40007	SCREW TO ADJUST CABLE TENSION	\$2.50
P40008	BEARING FOR FRONT SWIVEL FORK	\$3.50
P40011	NUT FOR FORK STEM P40035	\$2.50
P40012	STEEL WASHER FOR BRAKE SHOE	\$2.50
P40013 (PAIR)	CRIMP FOR END OF STEEL CABLE	\$2.50
P40014	WHITE PLASTIC WASHER FOR BRAKE SHOE	\$2.50
P40015	WHEEL NUT COVER	\$2.50
P40016	NUT FOR FRONT AXLE	\$2.50
P40023	WASHER FOR REAR AXLE	\$2.50
P40026	WHITE PLASTIC CABLE GUIDE - PART HELD ON BY RIVET	\$3.50
P40028	STEEL NUT FOR FOLDING MECH	\$2.50
P40032	FOLDING MECHANISM	\$6.95
P40035	SCREW TO ATTACH FORK TO FRAME	\$6.95
P40045	CENTER SCREW AND WASHER FOR FOLDING MECHANISM	\$2.50
P40047	REFLECTOR FOR B14-C	\$2.50
P42007	SCREW FOR FOLDING MECHANISM	\$2.50
P42009	CONTOURED WASHER FOR FOLDING MECHANISM	\$2.50
P42010	PLASTIC CONTOURED WASHER FOR FOLDING MECHANISM	\$2.50
P42012	NUT FOR BRAKE SHOE SCREW SHOE (SCREW LOCKS CABLE	\$2.50
	AND HOLDS SPRING)	
P42013 (PAIR)	BRAKE SPRING	\$3.50
P42014 (PAIR)	BRAKE SHOE	\$9.95
P42015	SCREW FOR BRAKE SHOE SHOE (SCREW LOCKS CABLE	\$2.50
	AND HOLDS SPRING)	
P42016	PLASTIC CONTOURED WASHER FOR REAR FORK	\$2.50
P42021	FORK FOR 6" FRONT WHEEL W/STEM	\$19.95
P42046	AXLE FOR FRONT WHEEL	\$2.50
P43003	LEFT SIDE RUBBER HAND GRIP	\$3.50
P43004	RIGHT SIDE RUBBER HAND GRIP	\$3.50
P-CABLE (PAIR)	CABLE FOR HAND BRAKE	\$3.50
P-COVER	BLACK CABLE COVER	\$3.50
V42007	U-SHAPE SEAT ATTACHMENT	\$2.50



V42009	A-FRAME SCREW SET FOR FRAME	\$2.50
V42052D	6" FRONT AND REAR WHEEL	\$19.95
V42064	SCREW & NUT FOR BACKREST HINGE	\$2.50
V42065	BLUE FLIP UP BACKREST FOR 4203BL	\$19.95
V42066	RED FLIP UP BACKREST FOR 4203RD	\$19.95
V42067	GREEN FLIP UP BACKREST FOR 4203GN	\$19.95
V42073	PLASTIC SEAT	\$19.95
V42073H	ANTI RATTLE CLAMP FOR PLASTIC SEAT	\$3.50
	(SCREWS INCLUDED)	
V42076	BLACK PLASTIC HINGE FOR BACKREST	\$2.50
V42079	LEFT HANDLE BAR	\$9.95
V42080	RIGHT HANDLE BAR	\$9.95
V42089	TRIANGLE SHAPED HEIGHT ADJUSTMENT KNOB	\$6.95
V42093	FOAM CUSHION FOR FLIP UP BACKREST	\$6.95
V42094	SLANTED A-FRAME PLUG	\$6.95
V42097	AXLE FOR REAR WHEEL	\$2.50
V42098	PINK FLIP UP BACKREST FOR 4203PK	\$19.95
V42099	BLACK FLIP UP BACKREST FOR 4203BK	\$19.95

**EXAMPLE SERIAL NUMBER 1106040203V4206AA0004551
AND SERIAL NUMBERS BEGINNING WITH "NV"**

NOVA Customer Service: 800-557-6682
www.NovaMedicalProducts.com