

**PLEASE READ BEFORE SEARCHING FOR THE PART NEEDED**



Many of the parts pictured on this breakdown may look identical to the part or parts needed but the measurements may be different and will not be interchangeable. Please check the serial number on your transport chair and make sure it corresponds with the serial number on this part diagram highlighted in yellow. The serial numbers on our transport chairs are located on a sticker on the crossbar under the seat as shown in the picture below. If you cannot find the serial number, please call customer service for assistance at (800) 557-6682 option 1.

Thank you



**NOVA Customer Service: 800-557-6682**

[www.NovaMedicalProducts.com](http://www.NovaMedicalProducts.com)

ITEM NUMBER 332 EXTRA WIDE LIGHTWEIGHT ALUMINUM TRANSPORT CHAIR W/HAND BRAKES

EXAMPLE SERIAL NUMBERS 1104352B0037 / 1104352R0037  
OR SERIAL NUMBER BEGINNING WITH "NV"



B14-C



CH-31929 CK-31903



CT-SB100



P-CABLE



P40007 P40013 PVT-33215



PVT-AC101



PVT-AN100 PVT-AN102



PVT-BS100



PVT-BS101



PVT-CS101



PVT-FC100



PVT-FW100



PVT-FW100A



PVT-HS100



PVT-PL100



PVT-PL101



PVT-RL101 PVT-WA100 PVT-WA101 PVT-WB102



S5000AH



V-33004



V33001



V33002



V33002A



V33003



V42010



VT-33201



VT-33202



VT-33203



VT-33204



VT-33205



VT-33206



VT-33207



VT-33208



VT-AP100



VT-AS100



VT-AS101



VT-BL100



VT-BR100



**VT-COVER**



**VT-F100B**



**VT-F100R**



**VT-FF100**



**VT-FP101**



**VT-FP102**



**VT-HS101**



**VT-P100**



**VT-SG100**



**VT-SU100**





VT-SU101



VT-TK101



VT-UB102



VT-US100



VT-WS100



**NOVA Customer Service: 800-557-6682**  
[www.NovaMedicalProducts.com](http://www.NovaMedicalProducts.com)

|                         |  |   |
|-------------------------|--|---|
| <b>ITEM NUMBER 332</b>  | <b>EXTRA WIDE LIGHTWEIGHT ALUMINUM TRANSPORT CHAIR W/HAND BRAKES</b> |  |
|                         | <b>EXAMPLE SERIAL NUMBERS 1104352B0037 / 1104352R0037</b>            |   |
|                         | <b>OR SERIAL NUMBER BEGINNING WITH "NV"</b>                          |   |
| <b>ITEM/PART NUMBER</b> | <b>DESCRIPTION</b>   |   |
| B14-C                   | HAND BRAKE WITH CABLE  |   |
| CH-31929                | BEARING & SPACER FOR FRONT WHEEL                                     | SET   |
| CK-31903                | TIP FOR BRAKE  |   |
| CT-SB100                | SEAT BELT  |   |
| P-CABLE                 | BRAKE CABLE  | PAIR  |
| P40007                  | TENSION SCREW FOR HAND BRAKE   |   |
| P40013                  | CRIMP FOR BRAKE CABLE  |   |
| PVT-33215               | ARM PAD (INCLUDES SCREWS)  |   |
| PVT-AC101               | COVER FOR REAR WHEEL   |   |
| PVT-AN100               | NUT FOR FRONT WHEEL AXLE   |   |
| PVT-AN102               | NUT FOR REAR WHEEL AXLE  |   |
| PVT-BS100               | SCREW SET FOR WHEEL BRAKE  | SET   |
| PVT-BS101               | SCREWS FOR B14-C HAND BRAKE  | PAIR  |
| PVT-CS101               | SCREW FOR CROSSBAR (INCLUDES WASHERS, NUT & SPACERS)                 | SET   |
| PVT-FC100               | COVER FOR FRONT SWIVEL FORK NUT/BEARINGS                             |   |
| PVT-FW100               | FRONT 8" WHEEL BLACK RUBBER TIRE (INCLUDES BEARINGS)                 |   |
| PVT-FW100A              | FRONT 8" WHEEL GRAY RUBBER TIRE (INCLUDES BEARINGS)                  |   |
| PVT-HS100               | SCREW FOR FOOT PLATE/HEEL STRAP (INCLUDES NUTS)                      |   |
| PVT-PL100               | PLUG FOR SEAT BAR  | PAIR  |
| PVT-PL101               | PLUG FOR FRAME   | PAIR  |
| PVT-RL101               | RELEASE LEVER FOR FLIP BACK ARM (INCLUDES SCREW & NUT)               |   |
| PVT-WA100               | AXLE FOR FRONT WHEEL   |   |
| PVT-WA101               | AXLE FOR REAR WHEEL  |   |
| PVT-WB102               | BEARING AND SPACER FOR 12" REAR WHEEL                                | SET   |
| S5000AH                 | RECEIVER FOR FLIP BACK ARM (INCLUDES SCREW & NUT)                    | PAIR  |
| V-33004                 | SEAT CRADLE  |   |
| V33001                  | KNOB FOR FOOT REST   |   |
| V33002                  | REAR 12" WHEEL BLACK RUBBER TIRE (INCLUDES BEARINGS)                 |   |
| V33002A                 | REAR 12" WHEEL GRAY RUBBER TIRE (INCLUDES BEARINGS)                  |   |
| V33003                  | SCREW FOR ARM PAD  | PAIR  |
| V42010                  | HAND GRIP  | PAIR  |
| VT-33201                | FLIP BACK DESK ARM LEFT IN BLUE                                      |   |
| VT-33202                | FLIP BACK DESK ARM RIGHT IN BLUE                                     |   |
| VT-33203                | FLIP BACK DESK ARM LEFT IN RED                                       |   |
| VT-33204                | FLIP BACK DESK ARM RIGHT IN RED                                      |   |
| VT-33205                | SCREW AND WASHER FOR SKIRT GUARD                                     |   |
| VT-33206                | SCREW, WASHER, AND NUT FOR SEAT CRADLE                               |   |



|           |  |      |
|-----------|--|------|
| VT-33207  | BEARING FOR FRONT SWIVEL FORK                        |      |
| VT-33208  | SPACER FOR FRONT SWIVEL FORK                         |      |
| VT-AP100  | SKIRT GUARD  |      |
| VT-AS100  | SUPPORT BRACKET FOR RIGHT FLIP BACK ARM              |      |
| VT-AS101  | SUPPORT BRACKET FOR LEFT FLIP BACK ARM               |      |
| VT-BL100  | BRAKE LEFT FOR REAR WHEEL                            | SET  |
| VT-BR100  | BRAKE RIGHT FOR REAR WHEEL                           | SET  |
| VT-COVER  | COVER FOR BRAKE CABLE                                |      |
| VT-F100B  | FOOT REST BLUE                                       | PAIR |
| VT-F100R  | FOOT REST RED  | PAIR |
| VT-FF100  | FORK FRONT SWIVEL                                    |      |
| VT-FP101  | FOOT PLATE RIGHT (COMPLETE)                          |      |
| VT-FP1012 | FOOT PLATE LEFT (COMPLETE)                           |      |
| VT-HS101  | HEEL STRAP   |      |
| VT-P100   | PIN FOR FOOT REST (INCLUDES PLASTIC WASHER AND NUT)  | PAIR |
| VT-P101   | RELEASE PIN FOR FOLDING MECHANISM                    |      |
| VT-SG100  | SEAT GLIDE   | PAIR |
| VT-SU100  | SCREWS AND WASHERS FOR SEAT UPHOLSTERY               | SET  |
| VT-SU101  | SCREWS AND WASHERS FOR BACK UPHOLSTERY               | SET  |
| VT-TK101  | TOE KICK   |      |
| VT-UB102  | BACK UPHOLSTERY (INCLUDES SCREWS AND WASHERS)        |      |
| VT-US102  | SEAT UPHOLSTERY (INCLUDES SCREWS, WASHERS, & STRAPS) |      |
| VT-WS100  | SPACER FOR REAR WHEEL AXLE                           |      |

**NOVA Customer Service: 800-557-6682**

[www.NovaMedicalProducts.com](http://www.NovaMedicalProducts.com)