

**PLEASE READ BEFORE SEARCHING FOR THE PART NEEDED**



Many of the parts pictured on this breakdown may look identical to the part or parts needed but the measurements may be different and will not be interchangeable. Please check the serial number on your transport chair and make sure it corresponds with the serial number on this part diagram highlighted in yellow. The serial numbers on our transport chairs are located on a sticker on the crossbar under the seat as shown in the picture below. If you cannot find the serial number, please call customer service for assistance at (800) 557-6682 option 1.

Thank you



**NOVA Customer Service: 800-557-6682**

[www.NovaMedicalProducts.com](http://www.NovaMedicalProducts.com)

ITEM NUMBER 349 LIGHTWEIGHT ALUMINUM TRANSPORT CHAIR W/FLIP ARMS

EXAMPLE SERIAL NUMBERS 1104352B0037 - 1104352R0037  
OR SERIAL NUMBER BEGINNING WITH "NV"



CK-31903



CT-SB100



P-CABLE



P40007



P40013



PVT-33215



PVT-AC102



PVT-AN100



PVT-AN102



PVT-CS102



PVT-FC100



PVT-FM100



PVT-FM101



**PVT-FW100A**



**PVT-HS100**



**PVT-PL100**



**PVT-PL101**



**PVT-RL100**



**PVT-RL101**



**PVT-RW100A**



**PVT-WA100 PVT-WA102**



**PVT-WB101**



S5000AH



V-33004



V33001



V33003



VT-33201



VT-33202



VT-33203



VT-33204



VT-33205



**VT-33206**



**VT-33207**



**VT-33208**



**VT-AP100**



**VT-AS100**



**VT-AS101**



**VT-BL101**



**VT-BR101**



**VT-COVER**



**VT-F100B**



**VT-F100R**



**VT-FF100**



**VT-FP101**



**VT-FP102**



**VT-HG100**



**VT-HS101**



**VT-LF106B**



**VT-LF106R**



**VT-P100**

**VT-P101**



**VT-RF106B**



**VT-RF106R**





VT-SG100



VT-SU100



VT-SU101



VT-TK101



VT-UB104




VT-US100



VT-WS100



**NOVA Customer Service: 800-557-6682**  
[www.NovaMedicalProducts.com](http://www.NovaMedicalProducts.com)

ITEM NUMBER 349	LIGHTWEIGHT ALUMINUM TRANSPORT CHAIR W/FLIP ARMS	
	<b>EXAMPLE SERIAL NUMBERS 1104352B0037 - 1104352R0037</b>	
	<b>OR SERIAL NUMBER BEGINNING WITH "NV"</b>	
ITEM/PART NUMBER	DESCRIPTION	
CK-31903	TIP FOR BRAKE	
CT-SB100	SEAT BELT	
P-CABLE	BRAKE CABLE	PAIR
P40007	TENSION SCREW FOR HAND BRAKE	
P40013	CRIMP FOR BRAKE CABLE	
PVT-33215	ARM PAD (INCLUDES SCREWS)	
PVT-AC102	COVER FOR REAR WHEEL	
PVT-AN100	NUT FOR FRONT WHEEL AXLE	
PVT-AN102	NUT FOR REAR WHEEL AXLE	
PVT-CS102	SCREW, SPACERS, AND WASHERS FOR CROSSBAR	SET
PVT-FC100	COVER FOR FRONT SWIVEL FORK NUT/BEARINGS	
PVT-FM100	FOLDING MECHANISM LEFT SIDE	
PVT-FM101	FOLDING MECHANISM RIGHT SIDE	
PVT-FW100A	FRONT 8" WHEEL GRAY RUBBER TIRE (INCLUDES BEARINGS)	
PVT-HS100	SCREW FOR FOOT PLATE/HEEL STRAP (INCLUDES NUTS)	
PVT-PL100	PLUG FOR SEAT BAR	PAIR
PVT-PL101	PLUG FOR FRAME	PAIR
PVT-RL100	LEVER FOR FOLDING MECHANISM	
PVT-RL101	RELEASE LEVER FOR FLIP BACK ARM (INCLUDES SCREW & NUT)	
PVT-RW100A	REAR 8" WHEEL GRAY RUBBER TIRE (INCLUDES BEARINGS)	
PVT-WA100	AXLE FOR FRONT WHEEL	
PVT-WA101	AXLE FOR REAR WHEEL	
PVT-WB101	BEARING FOR FRONT AND REAR WHEEL	
S5000AH	RECEIVER FOR FLIP BACK ARM (INCLUDES SCREW & NUT)	PAIR
V-33004	SEAT CRADLE	
V33001	KNOB FOR FOOT REST	
V33003	SCREW FOR ARM PAD	PAIR
VT-33201	FLIP BACK DESK ARM LEFT IN BLUE	
VT-33202	FLIP BACK DESK ARM RIGHT IN BLUE	
VT-33203	FLIP BACK DESK ARM LEFT IN RED	
VT-33204	FLIP BACK DESK ARM RIGHT IN RED	
VT-33205	SCREW AND WASHER FOR SKIRT GUARD	
VT-33206	SCREW, WASHER, AND NUT FOR SEAT CRADLE	
VT-33207	BEARING FOR FRONT SWIVEL FORK	
VT-33208	SPACER FOR FRONT SWIVEL FORK	
VT-AP100	SKIRT GUARD	
VT-AS100	SUPPORT BRACKET FOR RIGHT FLIP BACK ARM	



VT-AS101	SUPPORT BRACKET FOR LEFT FLIP BACK ARM	
VT-BL101	BRAKE LEFT FOR REAR WHEEL	SET
VT-BR101	BRAKE RIGHT FOR REAR WHEEL	SET
VT-COVER	COVER FOR BRAKE CABLE	
VT-F100B	FOOT REST BLUE	PAIR
VT-F100R	FOOT REST RED	PAIR
VT-FF100	FORK FRONT SWIVEL	
VT-FP101	FOOT PLATE RIGHT (COMPLETE)	
VT-FP1012	FOOT PLATE LEFT (COMPLETE)	
VT-HG100	GRIPS	PAIR
VT-HS101	HEEL STRAP	
VT-LF106B	HANDLE LEFT BLUE (INCLUDES FOLDING MECHANISM & GRIP)	SET
VT-LF106R	HANDLE LEFT RED (INCLUDES FOLDING MECHANISM & GRIP)	SET
VT-P100	PIN FOR FOOT REST (INCLUDES PLASTIC WASHER AND NUT)	PAIR
VT-P101	RELEASE PIN FOR FOLDING MECHANISM	
VT-RF106B	HANDLE RIGHT BLUE (INCLUDES FOLDING MECHANISM & GRIP)	SET
VT-RF106R	HANDLE RIGHT RED (INCLUDES FOLDING MECHANISM & GRIP)	SET
VT-SG100	SEAT GLIDE	PAIR
VT-SU100	SCREWS AND WASHERS FOR SEAT UPHOLSTERY	SET
VT-SU101	SCREWS AND WASHERS FOR BACK UPHOLSTERY	SET
VT-TK101	TOE KICK	
VT-UB104	BACK UPHOLSTERY (INCLUDES SCREWS AND WASHERS)	SET
VT-US100	SEAT UPHOLSTERY (INCLUDES SCREWS, WASHERS, & STRAPS)	SET
VT-WS100	SPACER FOR REAR WHEEL AXLE	

**NOVA Customer Service: 800-557-6682**

[www.NovaMedicalProducts.com](http://www.NovaMedicalProducts.com)