

PLEASE READ BEFORE SEARCHING FOR THE PART NEEDED



Many of the parts pictured on this breakdown may look identical to the part or parts needed but the measurements may be different. Please check the serial number on your rolling walker and make sure it corresponds with the serial number on this part diagram highlighted in yellow. The serial number on the 4010 rolling walker is located on a sticker under the left side crossbar as shown in the picture below. If you cannot find the serial number, please call customer service for assistance at (800) 557-6682 option 1.

Thank you




ITEM NUMBER 4010 CRUISER DE-LIGHT FOLDING WALKER

EXAMPLE SERIAL NUMBER 20110329003900100029

OR SERIAL NUMBER BEGINNING WITH "NT"



4010B



B15-C



CH-902608



P-CABLE



P-COVER



P40006



P40007



P40013



P40016



P40026



P40035



P40047



P40100



P40101



P40102



P40103



P40103B



P40103R



P40104



P40104B



P40104R



P40105LB

P40105LR

P40105RB

P40105RR



P40106L

P40106R

P40107

P40109

P40110



P42019

P42022

P42046



P42058



P401010

P401011

P401012



P401013




VT-33207



CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

ITEM NUMBER 4010	CRUISER DE-LIGHT FOLDING WALKER	
	EXAMPLE SERIAL NUMBER 20110329003900100029	
	OR SERIAL NUMBER BEGINNING WITH "NT"	
ITEM/PART NUMBER	DESCRIPTION	
4010B	BASKET	
B15-C	HAND BRAKE WITH CABLE	SET
CH-902608	PUSH BUTTON AND SPRING FOR LEGS	
P-CABLE	BRAKE CABLE	PAIR
P-COVER	COVER FOR BRAKE CABLE	
P40006	BEARINGS FOR WHEELS	SET OF 5
P40007	TENSION SCREW FOR HAND BRAKE	
P40013	CRIMP FOR BRAKE CABLE	PAIR
P40016	NUT FOR FRONT/REAR WHEEL AXLE	SET OF 5
P40026	CABLE GUIDE	
P40035	FORK STEM	
P40047	REFLECTOR FOR BRAKE HANDLE	
P40100	BARREL SCREW	
P40101	SCREW AND NUT FOR BRAKE SHOE	
P40102	SEAT	
P40103	SCREW AND NUT FOR HAND BRAKE	PAIR
P40103B	LEG REAR LEFT BLUE (COMPLETE WITH BRAKE HARDWARE AND WHEEL)	SET
P40103R	LEG REAR LEFT RED (COMPLETE WITH BRAKE HARDWARE AND WHEEL)	SET
P40104	FRAME CONNECTOR FOR CROSSBAR (INCLUDES HARDWARE)	SET
P40104B	LEG REAR RIGHT BLUE (COMPLETE WITH BRAKE HARDWARE AND WHEEL)	SET
P40104R	LEG REAR RIGHT RED (COMPLETE WITH BRAKE HARDWARE AND WHEEL)	SET
P40105LB	LEG FRONT WITH FORK AND WHEEL LEFT BLUE	SET
P40105LR	LEG FRONT WITH FORK AND WHEEL LEFT RED	SET
P40105RB	LEG FRONT WITH FORK AND WHEEL RIGHT BLUE	SET
P40105RR	LEG FRONT WITH FORK AND WHEEL RIGHT RED	SET
P40106L	BRAKE SHOE LEFT	
P40106R	BRAKE SHOE RIGHT	
P40107	AXLE FOR REAR WHEEL	
P40109	SPRING FOR BRAKE SHOE	
P40110	SCREW, NUT AND WASHER FOR SEAT	SET
P42019	WHEEL	
P42022	FORK FRONT SWIVEL (INCLUDES BEARING COVER AND BUSHING)	SET
P42046	AXLE FOR FRONT WHEEL	
P42058	PAD FOR SEAT	
P401010	CROSSBAR BLUE (INCLUDES CABLE GUIDE AND RIVET)	SET
P401011	KNOB FOLDING RELEASE	



P401012	CROSSBAR RED (INCLUDES CABLE GUIDE AND RIVET)	SET
P401013	COVER FOR BRAKE SHOE	
VT-33207	BEARING FOR SWIVEL FORK	

CONTACT US!
We're here to help



info@novajoy.com



1 800 557 6682