

PLEASE READ BEFORE SEARCHING FOR THE PART NEEDED



Many of the parts pictured on this breakdown may look identical to the part or parts needed but the measurements may be different. Please check the serial number on your rolling walker and make sure it corresponds with the serial number on this part diagram highlighted in yellow. The serial number on the 4239 and 4240 rolling walkers are located on a sticker under the seat on the crossbar the seat is attached to. Refer to the picture below. If you cannot find the serial number, please call customer service for assistance at (800) 557-6682 option 1.

Thank you



ITEM NUMBER 4239/4240 VIBE PETITE WIDE / VIBE WIDE

SERIAL NUMBER BEGINNING WITH "JN"



C40006



JN-420439



JN-423602



JN-423610



JN-423612



JN-423614



JN-423615



JN-423616



JN-423617



JN-423618



JN-423620



JN-423621



JN-423622



JN-423623



JN-423624



JN-423625



JN-423626



JN-423628



JN-423629



JN-423631



JN-423632



JN-423633



JN-423634



JN-423635



JN-423636



JN-423902



JN-423906



JN-423907



JN-423908



JN-423909



JN-423910



P40007



P40013



VT-33207



CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

ITEM NUMBER 4239/4240	VIBE PETITE WIDE / VIBE WIDE	
	SERIAL NUMBER BEGINNING WITH "JN"	
ITEM/PART NUMBER	DESCRIPTION	
C40006	SCREW TO ATTACH HAND BRAKE TO STEEL HANDLE	
JN-420439	REFLECTORS FOR HAND BRAKE HANDLE	PAIR
JN-423602	WHEEL FRONT 8" (INCLUDES BEARINGS)	
JN-423610	HINGE FOR BACKREST (INCLUDES HARDWARE)	SET
JN-423612	SEAT ATTACHMENTS (INCLUDES HARDWARE)	SET
JN-423614	KNOB WITH SCREW FOR HANDLE HEIGHT ADJUSTMENT	
JN-423615	HAND GRIPS	PAIR
JN-423616	HANDLE BAR LEFT	
JN-423617	HANDLE BAR RIGHT	
JN-423618	SCREWS FOR FOLDING MECHANISM (INCLUDES HARDWARE)	SET
JN-423620	BRAKES SHOES LEFT & RIGHT (INCLUDES HARDWARE)	SET
JN-423621	GUIDE FOR BRAKE CABLE	
JN-423622	HAND BRAKE WITH CABLE	SET
JN-423623	BRAKE CABLES	PAIR
JN-423624	COVER FOR BRAKE CABLES	
JN-423625	SPRING FOR BRAKE SHOE	
JN-423626	TENSION SCREW FOR BRAKE HANDLE	
JN-423628	AXLE FOR FRONT WHEEL (INCLUDES NUT)	SET
JN-423629	AXLE FOR REAR WHEEL (INCLUDES WASHERS & NUT)	SET
JN-423631	SILENCER FOR STEEL HANDLE BARS	
JN-423632	FORK FOR FRONT WHEEL (INCLUDES LOCK WASHER)	SET
JN-423633	BEARINGS & SPACERS FOR FRONT WHEEL	2 SETS
JN-423634	BEARINGS & SPACERS FOR REAR WHEEL	2 SETS
JN-423635	WHEEL REAR 8" (INCLUDES BEARINGS)	
JN-423636	PLUG SLANTED FOR A-FRAME	
JN-423902	BACKREST BLUE	
JN-423906	BACKREST RED	
JN-423907	CUSHION FOAM FOR BACKREST	
JN-423908	SEAT PADDED	
JN-423909	POUCH	
JN-423910	FOLDING MECHANISM	
P40007	TENSION SCREW FOR REAR LEG	
P40013	CRIMPS FOR BRAKE CABLE	PAIR
VT-33207	BEARING SWIVEL FOR FRONT FORK	

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682